**PROSPERIDAD
PARA TODOS**

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ANEXO 1

S0251 - Normal Priority Medical Device Alert

Medical Device

Special Report

Updated: January 9, 2014

UMDNS Terms:

- Infusion Pumps,
Multitherapy, Large
Volume, Multichannel
[17634]

Suggested Distribution:

- Clinical/Biomedical
Engineering
- IV Therapy
- Nursing

Geographic Regions:

- Worldwide

ECRI Institute User Experience Network—Lower Door Hinge May Crack on Carefusion Alaris 8100 Pump Modules

Product Identifier:

Model 8100 Alaris LVP Infusion Pump Modules [Capital Equipment]

Manufacturer:

- CareFusion Alaris 3750 Torrey View Ct, San Diego, CA 92130-2622,

Problem:

A member facility reports finding multiple CareFusion Alaris large-volume pump (LVP) modules with the lower left portion of the door cracked and chipped (see Figure 1). There are no reported injuries resulting from this problem.

Figure 1. Alaris Pump Module Door Hinge Chip

Discussion

CareFusion is currently investigating this problem. The damage occurs above the door's lower hinge. In most cases, the damage does not prevent the door from closing or latching. However, CareFusion has received a limited number of reports of the hinge breaking off when impacted. If the hinge breaks off, CareFusion recommends removing the module from service. In this scenario, a clinician may try to position the door so that it will latch; however, the lower portion of the door's inner platen may not be secure and flow accuracy may be compromised.

The member reports that 40% of their modules are affected but that no hinges have broken off. CareFusion provides replacement doors free of charge; the doors can be installed by clinical engineering. CareFusion indicates that this problem has been corrected in newly manufactured modules.

ECRI Institute Recommendations:

ECRI Institute recommends the following for users of the CareFusion Alaris System LVP:

- Do not use a LVP module in which the hinge has broken off. Remove the module from service until the door is replaced.
- Schedule repair for modules with cracked or chipped doors.

Comment:

- This alert is a living document and may be updated when ECRI Institute receives additional information. In circumstances in which we determine that it is appropriate for customers to repeat their review of an issue (e.g., when additional affected product has been identified), we will post a separate update alert. In other cases, we may add information, such as additional commentary, recommendations, and/or source documents, to the original alert. For additional information regarding the format of this alert, refer to our [HDA Format Guide](#).

Source(s):

- 2014 Jan 9. ECRI Institute researched member report.