ALERTA

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ANEXO

www.ecri.org . Printed from *Health Devices Alerts* on Friday, January 19, 2018 Page 1

[Critical Priority] - A29772 : Physio-Control—LIFEPAK 20e Defibrillator/Monitors: May Exhibit Power-Related Failures Medical Device Ongoing Action

Published: Monday, December 11, 2017
Last Updated: Thursday, December 14, 2017

UMDNS Terms:

- Defibrillator/Pacemakers, External [17882]

Product Identifier:

LIFEPAK 20e Defibrillator/Monitors [*Capital Equipment*]

For affected serial numbers, see the letter sent to your facility.

Units manufactured between September 2016 and June 2017

Geographic Regions: Worldwide

Manufacturer(s): Physio-Control Inc11811 Willows Rd NE, Redmond, WA 98052, United States

Suggested Distribution: Cardiology/Cardiac Catheterization Laboratory, Clinical/Biomedical Engineering, Critical Care, Emergency/Outpatient Services, OR/Surgery, Risk Management/Continuous Quality Improvement, Home Care, EMS/Transport

Problem:

In a December 2017 Urgent Medical Device Correction letter submitted by an ECRI Institute member hospital, Physio-Control states that it has received reports of the above devices exhibiting power-related failures as the user prepared the device for initial deployment or during use within the first year of distribution. The failures may include unexpected power on and power off, device lockup, or failure to power on or off. Any of these failures could result in a failure to deliver therapy to the patient and serious injury or death. Physio-Control also states that the failures are the result of a manufacturing process residue located beneath a component mounted on the power printed circuit board assembly (PCBA). Physio-Control has received 37 reports related to this problem; however, the firm has received no reports of adverse events related to the problem.

Action Needed:

Identify any affected product in your inventory. If you have affected product, verify that you have received the December 2017 Urgent Medical Device Correction letter and Confirmation Sheet from Physio-Control. A Physio-Control representative will contact your facility to arrange for a device correction that will include replacement of the power PCBA. Complete the Confirmation Sheet, and return it to Physio-Control using the instructions on the form.

For Further Information:

Physio-Control
Tel.: (866) 231-1220, 6 a.m. to 4 p.m. Pacific time, Monday through Friday
Website: [Click here](#)

Comments:

- This alert is a living document and may be updated when ECRI Institute receives additional information. In circumstances in which we determine that it is appropriate for customers to repeat their review of an issue (e.g., when additional affected product has been identified), we will post a separate update alert. In other cases, we may add information, such as additional commentary, recommendations, and/or source documents, to the original alert. For additional information regarding the format of this alert, refer to our [HDA Format Guide](#).

Source(s):

- 2017 Dec 8. Member Hospital. Physio-Control letter submitted by ECRI Institute member hospital [Download](#)
- 2017 Dec 14. Manufacturer. Manufacturer confirmed information