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ANEXO

www.ecri.org . Printed from *Health Devices Alerts* on Wednesday, September 6, 2017 Page 1

[High Priority] - A29148 : Siemens—ACUSON S Ultrasound Systems: Manufacturer Responds to WannaCry Ransomware Vulnerabilities Medical Device Ongoing Action

Published: Thursday, August 31, 2017

UMDNS Terms:

- Scanning Systems, Ultrasonic [14278]

Product Identifier:

ACUSON S Ultrasound Systems [Capital Equipment]

Software Version VE10A and All Software Versions of the following Software Releases: VC30, VC31, VD10

Geographic Regions: Impact in specific regions has not been identified or ruled out at the time of this posting), Worldwide

Manufacturer(s): Siemens Healthcare 40 Liberty Blvd, Malvern, PA, 19355, United States

Suggested Distribution: Cardiology/Cardiac Catheterization Laboratory, Clinical/Biomedical Engineering, Obstetrics/Gynecology/Labor and Delivery, Diagnostic Imaging, Information Technology

Problem:

In a Customer Advisory Notification letter submitted by an ECRI Institute member hospital, Siemens states that the Microsoft Windows operating system on the above systems has a vulnerability to the WannaCry virus. Siemens also states that if the system settings have been backed up and restored at any time, there is a high probability that the system is vulnerable to the virus. Siemens further states that exposure to the virus may halt system functionality, potentially resulting in a delay in procedure for a sedated patient or necessitating reinsertion of a transesophageal transducer or catheter, or repeat of a stress echo exam. Siemens states that patient harm as a result of this virus is improbable. For the above systems with software releases VC30 and VC31, a blue screen on the system indicates the likelihood that your network is infected. For the above systems with software release VD10 or version VE10A, Siemens has received no reports of infections; however, if the network becomes infected, there is a high probability that the system will also become infected. The manufacturer has not confirmed the information provided in the source material.

Action Needed:

Identify any affected systems in your inventory. If you have affected systems, verify that you have received the Customer Advisory Notification letter from Siemens. For systems with software releases VC30 and VC31, if you observe a blue screen, indicating that your network is infected, to continue to use the system you must disconnect the system from the infected network, restart the system, and then continue use of the system disconnected from the network. For affected systems with software release VD10 or version VE10A, if you observe symptoms of infection, such as lockups or ransomware messages, you must discontinue use of the system. Disconnect it from the network and contact your Siemens service representative. Siemens is initiating a product update to correct this vulnerability. For information on this update and others, click [here](#) and [here](#). Notify all relevant personnel at your facility of the information in the letter.

For further Information:

Siemens customer care center

Website: [Click here](#)

Comments:

- For information on other Siemens action related to the WannaCry ransomware, see Alerts [A28839](#), [A28842](#), and [A28925](#).
- For information on a Siemens action related to Microsoft Windows 7 and HP Client Automation vulnerabilities on various molecular imaging products, see [Alert A29058](#).
- This alert is a living document and may be updated when ECRI Institute receives additional information. In circumstances in which we determine that it is appropriate for customers to repeat their review of an issue (e.g., when additional affected product has been identified), we will post a separate update alert. In other cases, we may add information, such as additional commentary, recommendations, and/or source documents, to the original alert. For additional information regarding the format of this alert, refer to our [HDA Format Guide](#).