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ANEXO 1

www.ecri.org . Printed from *Health Devices Alerts* on Friday, July 10, 2015 Page 1

[Normal Priority] - S0276 : Draeger—Fabius GS Anesthesia Systems: Oxygen Flush Hose Cover Plate Required to Prevent Leaks [ECRI Exclusive User Experience Network] Medical Device Special Report

Published: Wednesday, June 24, 2015

UMDNS Terms:

- Anesthesia Units [10134]

Product Identifier:

Fabius GS Anesthesia Systems [Capital Equipment]

Geographic Regions: Worldwide

Manufacturer(s): Draeger Medical Inc 3135 Quarry Rd, Telford, PA 18969, United States

Suggested Distribution: Anesthesia, Clinical/Biomedical Engineering, OR/Surgery

Problem:

1. If the protective metal cover is not installed on the above systems (see Figure 1), the oxygen flush hose can deteriorate, which can cause the hose to leak. Units manufactured in 2006 or earlier are more likely to be missing the cover.
2. A significant leak in this hose would only affect the amount of oxygen delivered when pressing the oxygen flush button. Draeger states that a significant leak would be recognizable by the user during the manual part of the pre-use check. A non-significant leak in this hose will not affect the performance of the system.

Manufacturer's Corrective Action/Recommendations:

1. Draeger released a protective metal cover in 2006 that protects the hose and prevents deterioration. All units manufactured after 2006 should have the cover (see Figure 1).
2. A retrofit kit exists for units manufactured before 2006. Hospitals that are interested in the retrofit can contact their Draeger local representative to order the kit (part number MX08852; cost of the kit is approximately \$14). The kit can be installed by hospital staff; it does not require a Draeger service technician.

Figure 1. The oxygen flush hose on the Fabius GS without and with the cover.

ECRI Institute Recommendations:

Clinical/Biomedical Engineering

1. Determine if any Fabius GS units in your facility are missing the oxygen flush hose cover.
2. If any units do not have the cover, contact your Draeger local representative to order the retrofit kit and install the cover.

Comments:

- This alert is a living document and may be updated when ECRI Institute receives additional information. In circumstances in which we determine that it is appropriate for customers to repeat their review of an issue (e.g., when additional affected product has been identified), we will post a separate update alert. In other cases, we may add information, such as additional commentary, recommendations, and/or source documents, to the original alert. For additional information regarding the format of this alert, refer to our [HDA Format Guide](#).

Source(s):

- 2015 Jun 24. ECRI Institute researched member report.
- 2015 Jun 24. Figure 1. [Download](#)

©2015 ECRI Institute
5200 Butler Pike, Plymouth Meeting, PA 19462-1298, USA
May be reproduced by subscribing institution for internal distribution only.