## RECALL

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### ANEXO 1

www.ecri.org . Printed from *Health Devices Alerts* on Tuesday, December 13, 2016 Page 1

**[High Priority] - A27488 03 : \*Teleflex—LMA MAD Nasal Intranasal Mucosal Atomization Devices: May Not Deliver Fully Atomized Plume of Medication [Update]**  
Medical Device Ongoing Action

Published: Tuesday, November 29, 2016

#### UMDNS Terms:

- Atomizers [10226]

#### Product Identifier:

LMA MAD Nasal Intranasal Mucosal Atomization Devices [Consumable]

| Product Nos.: | Batch/Lot Nos.:                                                                                                          |
|---------------|--------------------------------------------------------------------------------------------------------------------------|
| MAD100        | 160105, 160137,<br>160302, 160321,<br>160402, 160435,<br>160506, 160523,<br>160609, 160620,<br>160707, 160802,<br>160813 |

|          |                           |
|----------|---------------------------|
| MAD100OS | 160322, 160524,<br>160630 |
|----------|---------------------------|

|        |                |
|--------|----------------|
| MAD110 | 160217, 160507 |
|--------|----------------|

|          |                |
|----------|----------------|
| MAD110OS | 160240, 160312 |
|----------|----------------|

|        |                           |
|--------|---------------------------|
| MAD130 | 160107, 160138,<br>160517 |
|--------|---------------------------|

|          |                |
|----------|----------------|
| MAD130OS | 160436, 160803 |
|----------|----------------|

|        |                                              |
|--------|----------------------------------------------|
| MAD140 | 160125, 160218,<br>160437, 160610,<br>160801 |
|--------|----------------------------------------------|

|          |                           |
|----------|---------------------------|
| MAD140OS | 160226, 160438,<br>160727 |
|----------|---------------------------|

|        |                                                                                                                                                                                                                                                                                                             |
|--------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| MAD300 | 160108, 160117,<br>160126, 160145,<br>160146, 160200,<br>160219, 160225,<br>160231, 160300,<br>160313, 160327,<br>160400, 160409,<br>160422, 160432,<br>160440, 160500,<br>160518, 160602,<br>160611, 160621,<br>160631, 160701,<br>160708, 160718,<br>160728, 160800,<br>160804, 160814,<br>160816, 160823 |
|--------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

©2016 ECRI Institute  
5200 Butler Pike, Plymouth Meeting, PA 19462-1298, USA  
May be reproduced by subscribing institution for internal distribution only.

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www.ecri.org . Printed from *Health Devices Alerts* on Tuesday, December 13, 2016 Page 2

MAD300B 160410

**Geographic Regions:** Argentina, Australia, Austria, Belgium, Canada, Chile, Colombia, Czech Republic, Denmark, &#160;Estonia, &#160;Finland, France, Germany, Hungary, Iceland, Iran, Ireland, Israel, &#160;Italy, Japan, Kuwait, &#160;Latvia, Luxembourg, Malta, The Netherlands, Philippines, Portugal, Reunion, Saudi Arabia, Slovenia, South Africa, Spain, Switzerland, United Arab Emirates, U.K., U.S.

**Manufacturer(s):** Teleflex Medical2917 Weck Drive Research Triangle Park, NC 27709, United States

**Suggested Distribution:** Critical Care, Emergency/Outpatient Services, Nursing, Pulmonology/Respiratory Therapy, Diabetes Education/Coordination, Neurology, Pain Clinic, EMS/Transport, Pharmacy, Materials Management

### Summary:

This Alert provides additional information based on a November 22, 2016, Amended Urgent Medical Device Recall Notification letter submitted by ECRI Institute member hospitals regarding [Alert Accession No. A27488](#) . Additional information is provided in the Problem and Action Needed fields.

### Problem:

[November 29, 2016]

In a November 22, 2016, Amended Urgent Medical Device Recall Notification letter submitted by ECRI Institute member hospitals, Teleflex states that in response to the October 27, 2016, Urgent Medical Device Recall Notification letter and the problem described below, some customers have indicated that, due to medical necessity, they plan to continue using the above devices rather than return them. Teleflex also states that it is advising all such customers to follow supplemental instructions for use (IFU) (attachment B of the [letter](#) ), which allow nondestructive testing of each unit before a procedure to determine if it is defective.

[October 28, 2016]

In an October 27, 2016, Urgent Medical Device Recall Notification letter submitted by an ECRI Institute member hospital, Teleflex states that the above devices may not deliver a fully atomized plume of medication. Teleflex also states that it has received reports of the above devices producing a straight stream instead of an atomized stream. Teleflex further states that failure of the device to deliver an atomized plume may impair the effectiveness of the medication with which it is used, potentially leading to serious injury or death in certain emergency situations, such as when the device is used in an off-label manner for needle-free delivery of drugs for reversal of life-threatening narcotic overdose, reversal of life-threatening hypoglycemia, or treatment of epileptic seizures.

### Action Needed:

Identify, isolate, and discontinue use of any affected devices in your inventory (if possible). If you have affected devices, verify that you have received the November 22, 2016, Amended Urgent Medical Device Recall Notification letter, supplemental IFU (attachment B), and Amended Recall Acknowledgment Form from Teleflex. Teleflex states that you should return affected devices for refund according to the instructions in the October 27, 2016, [letter](#) (for details, see the Action Needed field of [Alert Accession No. A27488](#) ). If you intend to continue using affected devices because of medical necessity, ensure that replacement devices are available and perform the following steps to nondestructively test each unit before use in a procedure, as instructed in the supplemental IFU attached to the [letter](#) :

- Attach a syringe containing 1 mL of either sterile water or sterile saline to the device.
- Briskly compress the plunger on the syringe so as to deliver the liquid through the device and observe how the liquid comes out at the (distal) end.
- If the liquid sprays in a fine mist then the device is atomizing as intended.
- If the device demonstrates streaming, select another MAD device for testing and use.

If your institution intends to continue using affected devices with the supplemental testing, complete the Amended Recall Acknowledgment Form and return it to Teleflex using the instructions on the form. Teleflex states that you may choose to follow either or both pathways (i.e., return for refund and/or continue to use with supplemental testing). For example, you may return lots 1 and 2 to Teleflex and notify the firm that you are retaining lot 3 and applying the supplemental instructions).

### For Further Information:

Teleflex local representative or Teleflex customer service department  
Tel.: (866) 246-6990

### Comments:

- For information on a similar and potentially related action initiated by Teleflex regarding LMA mucosal atomization devices, see [Alert Accession No. A27489](#) .
- This alert is a living document and may be updated when ECRI Institute receives additional information. In circumstances in which we determine that it is appropriate for customers to repeat their review of an issue (e.g., when additional affected product has been identified), we will post a separate update alert. In other cases, we may add information, such as additional commentary, recommendations, and/or source documents, to the original alert. For additional information regarding the format of this alert, refer to our [HDA Format Guide](#) .

### Source(s):

- 2016 Nov 28. Member Hospital. November 22, 2016, Teleflex letter submitted by an ECRI Institute member hospital (includes reply form) [Download](#)
- 2016 Nov 28. Member Hospital. October 27, 2016, Teleflex letter submitted by ECRI Institute member hospitals [Download](#)
- 2016 Nov 28. Manufacturer. Teleflex confirmed the information provided in the source material.