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ANEXO 1

www.ecri.org • Printed from Health Devices Alerts on Friday, May 30, 2014 Page 1

A22390 - High Priority Medical Device Alert

Medical Device Ongoing Action

Updated: May 29, 2014

UMDNS Terms:

- Ventilators, Intensive Care, Neonatal/Pediatric, High-Frequency [18793]

Suggested Distribution:

- Clinical/Biomedical Engineering
- NICU
- Pulmonology/Respiratory Therapy

Geographic Regions:

- Worldwide

CareFusion—3100A and 3100B High Frequency Oscillatory Ventilators: Manufacturer Issues Information Regarding Air Leak Patient Circuit Calibration and Performance Check

Product Identifier:

High Frequency Oscillatory Ventilators (HFOV): (1) 3100A, (2) 3100B [Capital Equipment]

Manufacturer:

- CareFusion Corp 3750 Torrey View Ct, San Diego, CA 92130, United States

Summary: On May 29, 2014, the manufacturer confirmed the information provided in the source material.

Problem: In an April 24, 2014, letter submitted by an ECRI Institute member hospital, CareFusion states that it has received reports of difficulty pressurizing patient circuits for the above ventilators during pre-use Patient Circuit Calibration. CareFusion also states that Patient Circuit Calibration and Performance Check are important to ensure the ventilator is ready for clinical use, and failure during these procedures may lead to delay in implementation of clinical treatment. Failure to perform pre-use checks as defined in the Operator's Manual may also cause the ventilator to pressurize incorrectly, potentially leading to damage to the ventilator and delay in treatment. For more information regarding troubleshooting Patient Circuit Calibration failures and common sources of circuit leaks, refer to the [letter](#) or the Operator's Manual (available [here](#)). CareFusion also states that contributing factors associated with a failure to pressurize, including the common sources of circuit leaks listed in the letter, may be misdiagnosed as the Cap Diaphragm leaking.

Action Needed:

Identify any affected product in your inventory. If you have affected product, verify that you have received the April 24, 2014, letter from CareFusion. CareFusion states that it is important that users refer to the Operator's Manual as needed before set-up and/or use of the above ventilators and that it is the user's responsibility to identify and eliminate all possible contributing factors associated with a failure to pressurize.

For Further Information:

CareFusion

Website: [Click here](#)

Comment:

- This alert is a living document and may be updated when ECRI Institute receives additional information. In circumstances in which we determine that it is appropriate for customers to repeat their review of an issue (e.g., when additional affected product has been identified), we will post a separate update alert. In other cases, we may add information, such as additional commentary, recommendations, and/or source documents, to the original alert. For additional information regarding the format of this alert, refer to our [HDA Format Guide](#).

Source(s):

- 2014 May 28. Member Hospital. CareFusion letter submitted by an ECRI Institute member hospital.
- 2014 May 29. Manufacturer. The manufacturer confirmed the information.